

Weekly Cleaning Procedure

CRYOTherm is constructed of the highest quality and most durable materials of any hydrotherapy device on the market. Proper care and maintenance of your CRYOTherm will extend the life of its components and provide a safe and sanitary environment for your athletes. The entire CRYOTherm should be completely drained, inspected, cleaned, and disinfected at least once a week using the following steps:

1. Turn all components of the control panel off.
2. Open the drain valve(s) to drain the tanks.
3. When tanks are empty remove the filter cartridge(s). See Changing the Filter Cartridge.
4. Remove any debris from around the pump suction, drain covers, filter housing, and hydrotherapy jets.
5. Use a hard surface disinfectant spray to clean all exposed surfaces of the tank, the Corian exterior, and the inside of the filter housing. Allow disinfectant to stand on the surfaces for at least 30 seconds before wiping.

Note: longer contact times may be required for certain pathogens, refer to product label for specific instructions.

6. Remove any stains in the tank - CRYOTherm has seamless stainless-steel tanks with coved corners which are easily wiped down and don't harbor bacteria. Stains that do not wipe away easily may be removed using Foaming Spray Cleaner and a Nylon Mesh Pad. Apply light to moderate pressure in the direction of the grain of the steel. Never use stainless steel cleaner to clean the tanks. Report any damage to your Grimm representative.
7. Remove stains on the Corian surface - the solid surface Corian exterior surface of CRYOTherm is also seamless, non-porous, and stain resistant. Any stubborn marks that do not easily wipe away can be buffed with a wet Nylon Mesh Pad, apply light to moderate pressure in a circular pattern.

Note: Only use Nylon Mesh Pads on white surfaces – do not use Nylon Mesh Pads on colored Corian surfaces such as lids and logos as it will cause scratches and swirling.

8. Install a clean filter cartridge(s), refer to Changing the Filter Cartridge.
9. Close drain valves.
10. Fill tanks according to procedure described in Fill Tanks.
11. Refer to Water Balance for information and advice about keeping water clean and sanitized.



Contact us at (800) 223-5395 or grimmscientific.com when you need assistance.