



Hydrotherapy,  
***Simplified***

**CRYO**Therm™

# CRYOTherm™

## ReadyFit

*Hydrotherapy solutions for every setting.*

The *ReadyFit* modular hydrotherapy system features popular tank sizes with dozens of configurations to fit your needs and your space.

This system is engineered to provide effective recovery for your athletes, and add phenomenal recruitment value to your athletic training room.

Use the [Build and Price Tool](#) on our website to build your optimal configuration.



*Every athlete deserves the ultimate recovery experience.*

# CustomFit

*Built for you and your space.*

When you need something more unique, *CustomFit* tubs are designed specifically for you. We can customize size, shape, and depth to fit any room.



	CustomFit™		
- Arrives Fully Operational	✓	✓	✓
- No Ice or Hot Water	✓	✓	✓
- No Daily Cleanings	✓	✓	✓
- Chemical Free	✓	✓	✓
- Cold or Cold <sub>2</sub> Warm	✓	✓	✓
- Custom Configuration	✓	✓	✓
- In-Ground Options		✓	✓
- FDA Listed Medical Device	✓	✓	✓
- Built to Local Health Dept Regulations			✓
- ETL Certification	✓	✓	
- Logo Inlay Option	✓	✓	✓

*Simple to install, easy to maintain,  
and built to last.*

## **PermitFit**

*When permits are required, PermitFit is the answer.*

We are able to engineer and build any chosen **CRYOTherm** configuration to be compliant with state and local health department and pool codes.

## **In-Ground**

*Uniquely designed for easier installation.*

Like all **CRYOTherm**s, our in-ground systems are fully assembled and tested in our factory, so they arrive on site ready for quick and easy installation. This lowers installation costs compared to other systems, making **CRYOTherm** the perfect choice for your in-ground project.




# CRYOTherm™

**CORIAN** surfaces are both durable and comfortable. They will provide years of maintenance free use.

CRYOTherm is always ready with **automatic water level** controls.

CRYOTherm features welded, seamless **stainless steel** tanks.

**COLD<sub>2</sub>WARM**  
Heat for the warm side is recovered from the cold side. Warm tank also switches to cold.



CORIAN inlaid logo available. **Logos** are a permanent part of the surface, not a sticker or decal.

CRYOTherm features a robust filtration and sanitation system to assure **clear and clean** water.

CRYOTherm has no standard tank size. **Customize** tank dimensions, seats, steps, and water depth.

**Precise settings** as low as 40°F for Cold tanks and as high as 110°F for Warm tanks.



*facility*  
**FRIENDLY**



*self*  
**CONTAINED**



*clinically*  
**EFFECTIVE**



*operationally*  
**EFFICIENT**

# FAQs

## **What does ‘Hydrotherapy, Simplified’ mean?**

We simplify the delivery of hydrotherapy, making your job more manageable. The **CRYOTherm** system is always full, always at the right temperature, always clean and sanitary, and always ready when your athletes need it. You drain and clean it once a week, and **CRYOTherm** does the rest, leaving you with more time for your other responsibilities.

## **How do I identify the best solution for my hydrotherapy needs?**

Call us! We’ve never found a hydrotherapy problem we couldn’t solve. We will conduct a site evaluation of your existing space and work with you to develop a plan. New construction? Tell us what you want and our engineering department will work with your architect directly to assure your hydrotherapy area looks and functions exactly the way you want it to.

## **How is my CRYOTherm installed?**

Installation supervision is offered with every **CRYOTherm** purchase. Our factory technician will deliver your hydrotherapy system and supervise its unloading and placement in the hydrotherapy area. Prior to production, we will ensure your unit fits your space and can be moved through doorways. The unit arrives completely self-contained so installation is quick and easy, even in difficult spaces. Our factory representative will also do a full in-service with your staff when the installation is complete.

## **What is the logo inlay?**

Logos are a permanent part of **CRYOTherm’s** exterior, not a sticker or decal. We will work with your branding department to develop an approved design consistent with your visual identity and hand craft the logo inlay using *Colors of Corian*. Logos create a stunning visual effect and are impactful for recruiting and facility appeal.

## Nationwide Service

Starting from day one and throughout the life of your **CRYOTherm**, Grimm is here to support you however we can. We provide online and phone support, as well as on-site service by one of our trained factory service personnel when necessary. It’s easy to get ahold of us and we do our best to meet your needs promptly.

