



Hydrotherapy,
Simplified

CRYOTherm™



*Every athlete deserves the
ultimate recovery experience.*



Exclusive with CRYOTherm

- **Stainless Steel Tanks**...ultimate durability and sanitation
- **CORIAN® Surfaces**...for long-lasting beauty and comfort
- **Robust Cooling System**...recovers heat to save energy
- **Water Level Controls**...shut off automatically when filling
- **Reliable Temperature Control**...accurate from 40°F to 104°F



facility
FRIENDLY



self
CONTAINED



clinically
EFFECTIVE




operationally
EFFICIENT


*Is hydrotherapy taking up
too much of your time?*

COLD₂WARM

CRYOTherm's proprietary heat recovery system utilizes energy removed from its cold tank to heat the warm tank. Warm tanks are also switchable to cold.



COLD



CRYOTherm's robust and reliable refrigeration system ensures stable cold plunge temperatures are maintained regardless of the number of athletes occupying its tank.

Our time-saving systems are always ready, without adding ice or hot water. You can trust **CRYO**Therm will be ready whenever your athletes need it.

Simple to install, easy to maintain and built to last.

Above-Ground

Least site prep—easiest to install and easily serviceable

CRYOTherms are designed to fit through standard doorways and, in most cases, use existing drain and utility connections.



Partial In-Ground

The look of In-Ground with the functionality of Above-Ground

Our partial in-ground systems are serviceable from the top so there is no climbing into a pit to change filters or drain tanks.



In-Ground

Uniquely designed for lower installation costs.

Like all **CRYO**Therms, our in-ground systems are fully assembled and tested in our factory, so they arrive on site ready for quick and easy installation.



What's Your Perfect Fit?

ReadyFit

Hydrotherapy solutions for every setting.

The **ReadyFit** modular hydrotherapy system features popular tank sizes with multiple configurations to fit your needs and your space.

CustomFit

Built for you and your space.

When you need something special, **CustomFit** systems are designed specifically for you. Customize size, shape, and features to fit any room.

PermitFit

When a permit is required, PermitFit is your answer.

Any chosen **CRYOTherm** configuration can be engineered for your state or local health department public pool code.

	<i>ReadyFit™</i>	<i>CustomFit™</i>	<i>PermitFit™</i>
Arrives Fully Operational	✓	✓	✓
No Ice or Hot Water	✓	✓	✓
No Daily Cleanings	✓	✓	✓
Built-In Sanitation	✓	✓	✓
Cold or Cold ₂ Warm	✓	✓	✓
Custom Configuration		✓	✓
In-Ground Options		✓	✓
FDA Listed Medical Device	✓	✓	✓
Built to Local Health Dept Regulations			✓
ETL Certification	✓	✓	✓
Logo Inlay Option	✓	✓	✓

Visit www.grimmscientific.com to identify models using the Build and Price Tool

We have the solution!

Easy to Install

Installation supervision is provided with every **CRYOTherm** purchase. Our factory technician will deliver your hydrotherapy system and supervise its unloading and placement in the hydrotherapy area.

Exclusive Logo Inlays

Corian logo inlays are a permanent part of **CRYOTherm**'s exterior, not a sticker or decal.

Hydrotherapy, Simplified


We simplify the delivery of hydrotherapy, making your job more manageable.

We are the best solution for your hydrotherapy needs

Call us! We've never found a hydrotherapy problem we couldn't solve.

Nationwide Service

Starting from day one and throughout the life of your **CRYOTherm**, Grimm is here to support you however we can. We provide online and phone support, as well as on-site service by one of our trained factory service personnel when necessary. It's easy to get ahold of us and we do our best to meet your needs promptly.



Visit www.grimmscientific.com or call (800) 223-5395 for more information.



Grimm Scientific Industries, Inc
PO Box 2143
Marietta, OH 45750

(800) 223-5395

www.grimmscientific.com

  @grimmscientific